VISCODERM[®]

Hydrogel
Patch

Trattamento lenitivo, idratante e rinfrescante
Soothing, hydrating and refreshing treatment

Contenuto:

1 maschera patch
per fronte, zigomi e area periorulare

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1 patch per collo/decolleté

Contains:

1 mask for forehead, cheek bones and perioral zone

1 mask for cheeks, chin and perioral zone

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INDICATIONS

VISCODERM® Hydrogel Patch formulation is indicated for soothing, refreshing and hydrating the skin after dermoesthetic treatments of the face, neck and décolleté.

Thanks to HYDROGEL technology, the ingredients are gradually released from the patch. The particular molecular structure of VISCODERM® Hydrogel Patch allows a significant quantity of water to be contained, which evaporates during the application, refreshing the treated area and releasing the ingredients in the skin through osmosis.

The new presentation of VISCODERM® Hydrogel Patch has been divided in 3, to easily and effectively cover each treated area.

INSTRUCTIONS FOR USE

Clean the area to be treated. Open the package and remove the specific patch required. Remove the PET protective film. Apply the patch making sure that it perfectly adheres to the skin surface without folds and preventing the formation of air bubbles. Keep patch applied on treated area for approximately 20 minutes in order to exert its action and after remove the patch. It's not necessary to rinse the face after the removal.

WARNINGS

For single-use only. Keep out of reach of children. Discontinue treatment in case of visible irritation or hypersensitivity.

The functionality and safety of the product are guaranteed only if the patch envelope is intact and undamaged. Once opened, the envelope must be resealed if all three patches are not used.

The product dries when exposed to air. External Use. Do not apply on damaged skin.

STORAGE

Keep away from direct light and heat sources. Dispose properly after use, do not abandon in the environment.

IBSA Farmaceutici Italia S.r.l.
Via Martiri di Cefalonia, 2
26900 Lodi – Italy
info@ibsaderma.com

MADE IN ITALY